



The all new **Thermasonic® Gel Warmer**

Add a little comfort and warmth to your patient's consultation experience.

Completely redesigned with microprocessor control, a powerful heating element, and higher sidewall heat exchangers to rapidly heat gel and consistently monitor temperature.



at the heart of each model

is advanced microprocessor technology and a Microtemp® thermal cutoff, providing protection from overheating. From basic, single bottle warming to degree-specific, three bottle temperature control, there is an option for almost any healthcare setting.

ThermoFisher
SCIENTIFIC

The world leader in serving science

THERMASONIC® GEL WARMER

NEW single bottle and two re-designed multi-bottle gel warmers rapidly warm ultrasound gel and consistently maintain temperature. Durable construction, suitable for wall mounting or counter top use and each unit comes with a two-year warranty. Options available include: Single Bottle with indicator, 3-Bottle holder with LED and 3-Bottle holder with LCD.

FEATURES

- Powerful heating element uniformly heats bottles of gel
- Microprocessor control maintains constant temperature
- "Thermal cutoff" provides protection from overheating
- Light-weight and portable design
- Low electrical demand
- High-density, impact-resistant composite
- Adjustable temperature settings for multi-bottle units
- Illuminated power ON/OFF switch on multi-bottle units
- Heat Indicator Lamp on single-bottle unit

Product #	Description	Pack Size
PKR82-01-20CE	Single-bottle gel warmer 230V	EA
PKR82-03-20CE	Multi-bottle gel warmer 230V with LEDs	EA
PKR83-03-20LCD	Multi-bottle gel warmer 230Vwith LCD display	EA



ThermoFisher
SCIENTIFIC

Australia: For customer service, call 1300-735-292
Visit us online at: www.thermofisher.com.au
New Zealand: For customer service, call 0800-933-966
Visit us online at: www.thermofisher.co.nz

©2015 Thermo Fisher Scientific Inc. All rights reserved.